



LiNA Electrodes

Single use and reusable electrodes for electrosurgery pencils



Accessories for pencils



E550-R-120
Electrode Extender.
For 2.35mm pencils.
Length 120mm

E550-R-200
Electrode Extender.
For 2.35mm pencils.
Length 200mm

Ordering Information:

LiNA Electrodes **0459**

REF Single use sterile	REF Reusable	Product	Width	Depth	Length	Packing
552-S-70	552-R-70	Blade electrode			70mm	Box of 10
552-S-70-C		Blade electrode coated			70mm	Box of 50
552-S-70-CI		Blade electrode coated & insulated			70mm	Box of 10
552-S-90	552-R-90	Blade electrode			90mm	Box of 10
552-S-125	552-R-125	Blade electrode			125mm	Box of 10
552-S-125-C		Blade electrode coated			125mm	Box of 10
552-S-125-CI		Blade electrode coated & insulated			125mm	Box of 10
552-S-165		Blade electrode			165mm	Box of 10
552-S-165-C		Blade electrode coated			165mm	Box of 10
552-S-165-CI		Blade electrode coated & insulated			165mm	Box of 10
552-S-200		Blade electrode			200mm	Box of 10
552-S-200-C		Blade electrode coated			200mm	Box of 10
552-S-200-CI		Blade electrode coated & insulated			200mm	Box of 10
553-S-70	553-R-70	Needle electrode			70mm	Box of 50
553-S-70-TU		Needle electrode, tungsten, 0,3mm			70mm	Box of 10
	554-R-3-120	Ball electrode. Diameter 3mm			120mm	Box of 10
554-S-5-120	554-R-5-120	Ball electrode. Diameter 5mm			120mm	Box of 10
C1315-S-120		Conization electrode	13mm	15mm	120mm	Box of 10
C1320-S-120		Conization electrode	13mm	20mm	120mm	Box of 10
C1608-S-120		Conization electrode	16mm	08mm	120mm	Box of 10
C1618-S-120		Conization electrode	16mm	18mm	120mm	Box of 10
C2020-S-120		Conization electrode	20mm	20mm	120mm	Box of 10
T1010-S-120	T1010-R-120	Loop T-bar electrode	10mm	10mm	120mm	Box of 10
T1512-S-120	T1512-R-120	Loop T-bar electrode	15mm	12mm	120mm	Box of 10
T2015-S-60		Loop T-bar electrode	20mm	15mm	60mm	Box of 10
T2015-S-120	T2015-R-120	Loop T-bar electrode	20mm	15mm	120mm	Box of 10
T2020-S-60		Loop T-bar electrode	20mm	15mm	60mm	Box of 10
T2020-S-120	T2020-R-120	Loop T-bar electrode	20mm	20mm	120mm	Box of 10
U2020-S-120		Loop U-bar electrode	20mm	20mm	120mm	Box of 10
	E550-R-120	Electrode Extender for 2.35mm pencils			120mm	Box of 10
	E550-R-200	Electrode Extender for 2.35mm pencils			200mm	Box of 10



Blade and Needle Electrodes



552-S-200 Blade electrode. Length 200mm - available in single use only



552-S-200-C Coated Blade electrode. Length 200mm - available in single use only



552-S-200-CI Coated/insulated Blade electrode. Length 200mm - available in single use only



552-S-165 Blade electrode. Length 165mm - available in single use only



552-S-165-C Coated Blade electrode. Length 165mm - available in single use only



552-S-165-CI Coated/insulated Blade electrode. Length 165mm - available in single use only



552-S-70 Blade electrode. Length 70mm



552-S-70-C Coated Blade electrode. Length 70mm - available in single use only



552-S-70-CI Coated/insulated Blade electrode. L 70mm - available in single use only



552-S-125 Blade electrode. Length 125mm



552-S-125-C Coated Blade electrode. Length 125mm - available in single use only



552-S-125-CI Coated/insulated Blade electrode. Length 125mm - available in single use only



553-S-70 Needle electrode. Length 70mm

Ball and Loop Electrodes

- For Cervical Biopsy, LLETZ & LEEP procedures



554-S-5-120 Ball electrode. Diameter 5mm. Length 120mm



T1010-S-120 T-bar. Width 10mm. Depth 10mm. Length 120mm



T1512-S-120 T-bar. Width 15mm. Depth 12mm. Length 120mm



T2015-S-120 T-bar. Width 20mm. Depth 15mm. Length 120mm



T2020-S-60 T-bar. Width 20mm. Depth 20mm. Length 60mm



T2020-S-120 T-bar. Width 20mm. Depth 20mm. Length 120mm



U2020-S-120 U-bar. Width 20mm. Depth 20mm. Length 120mm



T2015-S-60 T-bar. W 20mm. D 15mm. L 60mm

Conization Electrodes

- For Cervical Biopsy, LLETZ and LEEP procedures



C1618-S-60 Conization. W 16mm. D 18mm. L 60mm



C1608-S-120 Conization. Width 16mm. Depth 8mm. Length 120mm



C1315-S-120 Conization. Width 13mm. Depth 15mm. Length 120mm



C1618-S-120 Conization. Width 16mm. Depth 18mm. Length 120mm



C1320-S-120 Conization. Width 13mm. Depth 20mm. Length 120mm



C2020-S-120 Conization. Width 20mm. Depth 20mm. Length 120mm